

# Handheld Spirometer





# CHEST/MATE

CHESTMATE



**Quick,  
simple,  
accurate!!**

## Screen Shots



- Ideal for physician's office
- Portable and affordable
- No periodical calibration required
- Sophisticated and intuitive PC software
- Easy and fast data transmission

**CHEST**

CE 0123

## Specifications

Flow detection	Turbine flow sensor
Flow range	±18 L/s
Flow accuracy	±5% or ±0.2 L/s, whichever is greater
PEF accuracy	±10% or ±0.17 L/s, whichever is greater
Volume detection	Flow integration
Volume range	±10 L/s
Volume accuracy	±2.5% or ±0.05 L, whichever is greater
Display	3.5 inch color touchscreen monitor
Internal memory	10,000 data
Interface	USB (Micro USB-AB)
Print	Via PC software (Spirobank-02)
Internal battery	Lithium-ion
Dimensions	150x66x21 (mm)
Weight	125 g

Specifications are subject to change without notice.

## Optional Carrying Case



## Analyzed Items

FVC, FEV0.5, FEV0.75, FEV1.0, FEV2.0, FEV3.0, FEV6.0, FEV0.5/FVC, FEV0.75/FVC, FEV1.0/FEV6.0, FEV1.0/FVC, FEV2.0/FVC, FEV3.0/FVC, FEV6.0/FVC, FVC/HT, FEV1.0/HT, FET, ExtrapV, ExtrapV%, PEF, PEF TIME, FEF25-75%, FEF25%, FEF50%, FEF75%, FEF90%, FEF50%/FEF75%, PEF/HT, FEF25%/HT, FEF50%/HT, FEF75%/HT, FEF75-85%, FEF200-1200, MTC75-50, MTC50-25, MTC25-RV, MTCR, OI, ATI, MVV43, FVC+FEV1.0, FEF25-75%/FVC, FIVC, FIV0.5, FIV1.0, FIV1.0/FVC, FIV1.0/FIVC, PIF, FIF50%, FEF50/FIF50, FIF50/FEF50, CVI

## Standard Accessories

Turbine Flow Sensor Head	1pc
Paper Mouthpiece	1pc
Plastic Nose Clip	1pc
USB Cable (A-micro B)	1pc
Stylus Pen	1pc
PC Software (Spirobank-02)	Downloadable

## Optional Accessories


Carrying Case	1pc
Calibration Syringe (3L or 1L)	1pc
Paper Mouthpiece	100pcs/box
Spiro-Filter (Standard or Mouthpiece type)	100pcs/box
Plastic Nose Clip	10pcs/pack

Other optional accessories available.


## CHEST SpiroBank-02



All Results



Calibration Check



FVC Measurement



Trend Graph

## CHEST M.I.,INC.

3-25-11 Hongo Bunkyo-ku, Tokyo 113-0033, Japan  
 Phone : +81-3-5804-5031 Fax : +81-3-3812-7284  
<https://www.chest-mi.co.jp/en/>

DISTRIBUTED BY :